

2-Phenylbutyric acid


Catalog No: #S0631

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2-Phenylbutyric acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	164.2
Formulation	C10H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCC(C(O)=O)C1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. MATHIVAT A, COTTET J. Bull Mem Soc Med Hop Paris. 1953 Nov 27-Dec 4;69(30-31):1030-48.

Note: This product is for in vitro research use only