

## Gallic acid trimethyl ether

Catalog No: #S0632


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Gallic acid trimethyl ether
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	212.2
Formulation	C10H12O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: COC1=CC(=CC(OC)=C1OC)C(=O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Limited solubility  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Doutheil J, Paschen W. Pharmacol Toxicol. 1999 Aug;85(2):85-91.

Note: This product is for in vitro research use only