

(2-Chlorophenyl)phenyl-methanone


Catalog No: #S0633

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	(2-Chlorophenyl)phenyl-methanone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	216.66
Formulation	C13H9ClO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: Clc1ccccc1C(=O)c2ccccc2
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Sancho MI, et al. Carbohydr Res. 2011 Sep 27;346(13):1978-84.

Note: This product is for in vitro research use only