

## 2-Acetyl-4-butylamidophenol

Catalog No: #S0634


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2-Acetyl-4-butylamidophenol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	221.25
Formulation	C <sub>12</sub> H <sub>15</sub> NO <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CCCC(=O)NC1=CC(C(=O)C)=C(O)C=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Limited  
solubility  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Wiesmann UN, et al. Biochem Biophys Res Commun. 1975 Oct 27;66(4):1338-43.

Note: This product is for in vitro research use only