

## p-(2-Methoxyethyl) phenol


Catalog No: #S0635

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	p-(2-Methoxyethyl) phenol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	152.19
Formulation	C9H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
SMILES: COCCC1=CC=C(O)C=C1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: Soluble in Ethanol, slightly soluble in Water  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only