

2-Acetyl-6-methoxynaphthalene

Catalog No: #S0637


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-Acetyl-6-methoxynaphthalene
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	200.23
Formulation	C13H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: COC1=CC2=C(C=C1)C=C(C=C2)C(C)=O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt:
Solubility: Limited solubility
Appearance: White Needle Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kim J, et al. J Phys Chem B. 2015 Feb 12;119(6):2620-7.

Note: This product is for in vitro research use only