

## 2-Methyl-4-nitroimidazole

Catalog No: #S0638


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2-Methyl-4-nitroimidazole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	127.1
Formulation	C4H5N3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:ImmunologySMILES:CC1=NC=C(N1)[N+](O-)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Off-white Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Miller MW, et al. J Med Chem. 1970 Sep;13(5):849-52.

Note: This product is for in vitro research use only