

3-Amino-5-methylisoxazole


Catalog No: #S0644

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-Amino-5-methylisoxazole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	98.1
Formulation	C4H6N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC1=CC(N)=NO1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited solubility
Appearance: White
Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. M J Rieder et al. Journal of Pharmacology and Experimental Therapeutics, 244(2), 724-728 (1988)

Note: This product is for in vitro research use only