

Tryptamine

Catalog No: #S0647

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tryptamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor agonist; Dehydrogenase inhibitor; TRP/TRPV Channel inhibitor
Calculated MW	160.22
Formulation	C10H12N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: NCCC1=CNC2=CC=CC=C12
Pathways: GPCR/G Protein; Membrane transporter/Ion channel; Metabolism;
Neuroscience
Receptor: 5-HT; AADH; hTRPA1
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Faint Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Berman HM, et al. Nucleic Acids Res. 2000 Jan 1;28(1):235-42.

Note: This product is for in vitro research use only