

Phenyl acetate


Catalog No: #S0650

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phenyl acetate
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	136.15
Formulation	C8H8O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CC(=O)OC1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhang H, et al. PLoS One. 2014 May 5;9(5):e96518.

Note: This product is for in vitro research use only