

N-Phenylbenzylamine


Catalog No: #S0651

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-Phenylbenzylamine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Microtubule Associated inhibitor
Calculated MW	183.25
Formulation	C13H13N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C(NC1=CC=CC=C1)C1=CC=CC=C1
Pathways: Cytoskeletal Signaling
Receptor: Tubulin polymerization
Boiling pt: Melting pt: Solubility: Limited solubility
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Cushman M, et al. J Med Chem. 1993 Sep 17;36(19):2817-21.

Note: This product is for in vitro research use only