

Tris(2-chloroethyl) phosphate


Catalog No: #S0653

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tris(2-chloroethyl) phosphate
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	285.49
Formulation	C6H12Cl3O4P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1CCOP(=O)(OCCCl)OCCCl
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Water: 7,000 mg/L
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Cheng W, et al. Chemosphere. 2013 Aug;92(8):1017-21.

Note: This product is for in vitro research use only