

2,3-Dihydrofuran


Catalog No: #S0654

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2,3-Dihydrofuran
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	70.09
Formulation	C ₄ H ₆ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: C1C=CCO1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Morel A, et al. *Molecules*. 2014 Jun 19;19(6):8402-13.

Note: This product is for in vitro research use only