

3,4-Dimethoxybenzyl alcohol


Catalog No: #S0657

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3,4-Dimethoxybenzyl alcohol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	168.19
Formulation	C9H12O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=C(OC)C=C(CO)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tien M, et al. J Biol Chem. 1986 Feb 5;261(4):1687-93.

Note: This product is for in vitro research use only