

2-(2-Hydroxyethyl)pyridine


Catalog No: #S0658

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2-(2-Hydroxyethyl)pyridine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	123.15
Formulation	C7H9NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OCCC1=CC=CC=N1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: soluble in Ethanol
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Icel C, et al. Eur J Med Chem. 2013 Feb;60:386-94.

Note: This product is for in vitro research use only