

Bispyrazolone


Catalog No: #S0659

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bispyrazolone
Brief Description	Inhibitors
Purification	99.86%
Target Name	Others
Calculated MW	346.38
Formulation	C20H18N4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC1=NN(C(=O)C1C1C(C)=NN(C1=O)C1=CC=CC=C1)C1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only