

## Cyclin A1 Polyclonal Antibody

Catalog No: #40801

Package Size: #40801-1 50ul #40801-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cyclin A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Cyclin A1 Polyclonal Antibody detects endogenous levels of Cyclin A1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Cyclin A1.
Target Name	Cyclin A1
Other Names	CCNA1; Cyclin-A1
Accession No.	Swiss-Prot: P78396NCBI Gene ID: 8900
Uniprot	P78396
GeneID	8900;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of HeLa cells using Cyclin A1 Polyclonal Antibody



Western Blot analysis of Jurkat NIH-3T3 cells using Cyclin A1 Polyclonal Antibody

---

Note: This product is for in vitro research use only