

## Beta-Tetralone

Catalog No: #S0660


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Beta-Tetralone
Brief Description	Inhibitors
Purification	97.00%
Target Name	P450 antagonist
Calculated MW	146.19
Formulation	C10H10O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: O=C1CCC2=CC=CC=C2C1  
Pathways: Metabolism  
Receptor: CYP2A  
Boiling pt: Melting pt: Solubility: Soluble in DMSO  
Appearance: Liquid  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. ADKINS H, et al. J Am Chem Soc. 1948 Dec;70(12):4247.

Note: This product is for in vitro research use only