

Dimethylsulfamoyl chloride

Catalog No: #S0661


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dimethylsulfamoyl chloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	143.59
Formulation	C2H6ClNO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CN(C)S(Cl)(=O)=O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kevill DN, et al. Org Biomol Chem. 2006 Apr 21;4(8):1580-6.

Note: This product is for in vitro research use only