

Nicotinic acid hydrazide

Catalog No: #S0664


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nicotinic acid hydrazide
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	137.14
Formulation	C6H7N3O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: NNC(=O)C1=CN=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kruszynski R, et al. Acta Crystallogr C. 2011 Feb;67(Pt 2):o52-6.

Note: This product is for in vitro research use only