

2,6-Dihydroxypurine

Catalog No: #S0665

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2,6-Dihydroxypurine
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	152.11
Formulation	C ₅ H ₄ N ₄ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: O=C1NC2=C(NC(=O)N1)N=CN2
Pathways: Others
Receptor: PNP2; XPT
Boiling pt: Melting pt: Solubility: Water: 69 mg/L (at 16 C)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Overington JP, et al. Nat Rev Drug Discov. 2006 Dec;5(12):993-6.

Note: This product is for in vitro research use only