

4-Chlorosalicylic acid


Catalog No: #S0666

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Chlorosalicylic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	172.57
Formulation	C7H5ClO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC(=O)C1=C(O)C=C(Cl)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited
solubility Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wei Jiang et. al. Analytical Sciences, 19(6), 923-925 (2003)

Note: This product is for in vitro research use only