

Flecainide Impurity D

Catalog No: #S0667


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Flecainide Impurity D
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	318.17
Formulation	C11H8F6O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: OC(=O)C1=CC(OCC(F)(F)F)=CC(OCC(F)(F)F)=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting
 pt: Solubility: Limited solubility
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Cecile Terrenoire et al. Journal of General Physiology, 141(1), 61-72 (2013)

Note: This product is for in vitro research use only