

1-Aminohydantoin hydrochloride

Catalog No: #S0668


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	1-Aminohydantoin hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	151.55
Formulation	C3H6CIN3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:Cl.NN1CC(=O)NC1=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Limited solubilityAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Liu W, et al. J Agric Food Chem. 2007 Aug 22;55(17):6829-34.

Note: This product is for in vitro research use only