

## Chloramine-T

Catalog No: #S0669


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Chloramine-T
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	227.65
Formulation	C7H7ClNNaO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: [Na+].CC1=CC=C(C=C1)S(=O)(=O)[N-]Cl  
Pathways: Others  
Receptor: Others  
Boiling pt: 314.3°C  
Melting pt: 167-170°C  
Solubility: Soluble in Water  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Boran H, Altinok I. J Fish Dis. 2014 May; 37(5):431-41.

Note: This product is for in vitro research use only