

## L-arginine monohydrochloride

Catalog No: #S0670


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	L-arginine monohydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives antagonist
Calculated MW	210.66
Formulation	C6H15ClN4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
 SMILES: C(C[C@@H](C(=O)O)N)CN=C(N)N  
 Pathways: Others  
 Receptor: Amino Acids and Derivatives  
 Boiling pt: Melting pt: 226-230°C  
 Solubility: Water: 42 mg/mL (199.37 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Morris CR, et al. Am J Respir Crit Care Med, 2003, 168(1), 63-69.

Note: This product is for in vitro research use only