

2-Naphthol

Catalog No: #S0671


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-Naphthol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	144.17
Formulation	C10H8O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:OC1=CC2=CC=CC=C2C=C1Pathways:OthersReceptor:OthersBoiling pt:285-286CMelting pt:120-122CSolubility:DMSO: 10 mMAppearance:Grey SheetRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ad?o P, et al. Dalton Trans. 2015 Jan 28; 44(4):1612-26.

Note: This product is for in vitro research use only