

## Mevastatin

Catalog No: #S0683


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Mevastatin
Brief Description	Inhibitors
Purification	99.00%
Target Name	HMG-CoA Reductase inhibitor; Lipid inhibitor
Calculated MW	390.51
Formulation	C23H34O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: CC[C@H](C)C(=O)O[C@H]1CCC=C2[C@H]1[C@H]([C@H](C=C2)C)CC[C@@H]3C[C@H](CC(=O)O3)O Pathways: Metabolism Receptor: HMG-CoA reductase; CES1 Boiling pt: 55°C Melting pt: 151-153°C Solubility: DMSO: 78mg/ml (199.74 mM); Water:

## References

1. Endo A, et al. J Antibiot (Tokyo), 1976, 29(12), 1346-1348.

Note: This product is for in vitro research use only