

Retinyl acetate

Catalog No: #S0685


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinyl acetate
Brief Description	Inhibitors
Purification	95.00%
Target Name	Vitamin
Calculated MW	328.5
Formulation	C22H32O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CC1=C(C(CCC1)(C)C)/C=C/C(=C/C=C/C(=C/COC(=O)C)/C)/CPathways:MetabolismReceptor:VABoling pt:440.5CMelting pt:57-58CSolubility:Soluble in EthanolAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Stephensen CB, et al. Annu Rev Nutr, 2001, 21, 167-192.

Note: This product is for in vitro research use only