

## Ganciclovir

Catalog No: #S0688


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ganciclovir
Brief Description	Inhibitors
Purification	99.55%
Target Name	Antifection inhibitor
Calculated MW	255.23
Formulation	C <sub>9</sub> H <sub>13</sub> N <sub>5</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: NC1=NC(=O)C2=C(N1)N(COC(CO)CO)C=N2 Pathways: Microbiology&virology Receptor: FHV-1 Boiling pt: 675°C Melting pt: 250°C Solubility: DMSO: 50 mg/mL (195.9 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Champness JN, et al. Proteins. 1998 Aug 15;32(3):350-61.

Note: This product is for in vitro research use only