

Sulfamonomethoxine

Catalog No: #S0699


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfamonomethoxine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	280.31
Formulation	C ₁₁ H ₁₂ N ₄ O ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: COC1=NC=NC(NS(=O)(=O)C2=CC=C(N)C=C2)=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 204-206°C
Solubility: DMSO: 42mg/mL
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Ma J, et al. Environ Sci Pollut Res Int. 2014 Jun;21(12):7436-48.

Note: This product is for in vitro research use only