

## Broxiquinoline

Catalog No: #S0701


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Broxiquinoline
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	302.95
Formulation	C <sub>9</sub> H <sub>5</sub> Br <sub>2</sub> NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: OC1=C(Br)C=C(Br)C2=C1N=CC=C2  
Pathways: Others  
Receptor: Others  
Boiling pt: 370.5  
Melting pt: 198-200  
Solubility: DMSO: 17 mg/mL (56.11 mM)  
Appearance: Gray Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Fritzier JM, et al. J Antimicrob Chemother, 2012, 67(3), 609-617.

Note: This product is for in vitro research use only