

L-Cycloserine


Catalog No: #S0710

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-Cycloserine
Brief Description	Inhibitors
Purification	99.92%
Target Name	Others inhibitor
Calculated MW	102.09
Formulation	C3H6N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:C1[C@@H](C(=O)NO1)NPathways:OthersReceptor:OthersBoiling pt:Melting pt:147CSolubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Noda M, et al. J Biol Chem. 2004 Oct 29;279(44):46143-52.

Note: This product is for in vitro research use only