

Cyclin G Polyclonal Antibody

Catalog No: #40806

Package Size: #40806-1 50ul #40806-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin G Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Cyclin G Polyclonal Antibody detects endogenous levels of Cyclin G protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Cyclin G.
Target Name	Cyclin G
Other Names	CCNG1; CCNG; CYCG1; Cyclin-G1; Cyclin-G
Accession No.	Swiss-Prot: P51959NCBI Gene ID: 900
Uniprot	P51959
GeneID	900;
SDS-PAGE MW	29kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using Cyclin G Polyclonal Antibody

Note: This product is for in vitro research use only