

## Suloctidil

Catalog No: #S0712


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Suloctidil
Brief Description	Inhibitors
Purification	99.00%
Target Name	Platelet aggregation
Calculated MW	337.56
Formulation	C <sub>20</sub> H <sub>35</sub> NOS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: CCCCCCN[C@H](C)[C@H](C1=CC=C(C=C1)SC(C)C)O Pathways: Others Receptor: Antiplatelet-aggregation

Boling pt: 467.3°C Melting pt: Solubility: DMSO: 10 mM Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Ishibashi A, et al. Nihon Yakurigaku Zasshi. 1983 Nov;82(5):361-73.

Note: This product is for in vitro research use only