

Chlormezanone

Catalog No: #S0715


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chlormezanone
Brief Description	Inhibitors
Purification	99.89%
Target Name	Others
Calculated MW	273.73
Formulation	C ₁₁ H ₁₂ ClNO ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CN1C(C2=CC=C(Cl)C=C2)S(=O)(=O)CCC1=O
Pathways: Others
Receptor: Translocator protein
Boiling pt: 534.5°C
Melting pt: 114°C
Solubility: Ethanol: 26 mg/mL (94.98 mM); DMSO: 55 mg/mL (200.92 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Overington JP, et al. Nat Rev Drug Discov. 2006 Dec;5(12):993-6.

Note: This product is for in vitro research use only