

Gallamine triethiodide


Catalog No: #S0720

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gallamine triethiodide
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR antagonist
Calculated MW	891.53
Formulation	C30H60I3N3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CCN(CC)CCOC1=CC=CC(OCCN(CC)CC)=C1OCCN(CC)CC
Pathways: Neuroscience
Receptor: AChR
Boiling pt: 502.6°C
Melting pt: 235°C
Solubility: DMSO: 100mg/ml (112.17 mM); Water: 100mg/ml (112.17 mM); Ethanol:

References

1. Huang XP, et al. Mol Pharmacol. 2005 Sep;68(3):769-78.

Note: This product is for in vitro research use only