

## Piroctone olamine

Catalog No: #S0721


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Piroctone olamine
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	298.42
Formulation	C16H30N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine system SMILES:CC(CC1=CC(C)=CC(=O)N1O)CC(C)(C)C Pathways:Others Receptor:Others Boiling pt:344.1C? Melting pt:130-135C Solubility:DMSO: 10 mM Appearance:White Crystal Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Dubini F, et al. *Arzneimittelforschung*. 2005;55(8):478-83.

Note: This product is for in vitro research use only