

## Ursolic acid

Catalog No: #S0722


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ursolic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	456.7
Formulation	C30H48O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:CancerSMILES:C[C@@H]1CC[C@@]2(CC[C@@]3(C(=CC[C@H]4[C@@]3(CC[C@@H]5[C@@]4(CC[C@@H](C5(C)C)O)C)C)[C@@H]2[C@H]1C)C(=O)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:292CSolubility:DMSO: 91 mg/mL (199.25 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Kim M, et al. Int J Mol Med. 2015 Mar;35(3):755-62.

Note: This product is for in vitro research use only