

## Cloperastine hydrochloride

Catalog No: #S0723


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cloperastine hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	5-HT Receptor
Calculated MW	366.33
Formulation	C <sub>20</sub> H <sub>25</sub> Cl <sub>2</sub> NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Respiratory systemSMILES:Cl.C1C1=CC=C(C=C1)C(OCCN1CCCC1)C1=CC=CC=C1Pathways:GPCR/G Protein;  
NeuroscienceReceptor:HTBoling pt:423.9CMelting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility ,  
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Yu MC, et al. Eur J Pharmacol, 2010, 627(1-3), 269-275.

Note: This product is for in vitro research use only