

## Formononetin

Catalog No: #S0724


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Formononetin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	268.27
Formulation	C16H12O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:COC1=CC=C(C=C1)C1=COC2=C(C=CC(O)=C2)C1=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:256-260CSolubility:DMSO: 54 mg/mL (201.29 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Liu XJ, et al. Asian Pac J Y Prev. 2014;15(6):2835-9.

Note: This product is for in vitro research use only