

## Benzydamine hydrochloride

Catalog No: #S0726


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Benzydamine hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	345.87
Formulation	C19H23N3OHCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:Cl.CN(C)CCCOC1=NN(CC2=CC=CC=C2)C2=CC=CC=C12Pathways:Immunology/Inflammation;  
NeuroscienceReceptor:COXBoling pt:474.4C?Melting pt:147-153CSolubility:Ethanol: 69 mg/mL (199.49 mM); Water: 69 mg/mL (199.49 mM); DMSO:  
69 mg/mL (199.49 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the  
ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Epstein JB, et al. Int J Radiat Oncol Biol Phys, 1989, 16(6), 1571-1575.

Note: This product is for in vitro research use only