

Piperonyl butoxide

Catalog No: #S0727


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Piperonyl butoxide
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	338.44
Formulation	C19H30O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCCCOCCOCCOC1=CC=C(OCO2)C=C1CCCPathways:OthersReceptor:MFOsBoiling pt:Melting pt:180CSolubility:Soluble in WaterAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Casida, J. E. J Agric Food Chem. 1970 Sep-Oct;18(5):753-72.

Note: This product is for in vitro research use only