

Benzyl benzoate


Catalog No: #S0729

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Benzyl benzoate
Brief Description	Inhibitors
Purification	99.36%
Target Name	Others
Calculated MW	212.24
Formulation	C ₁₄ H ₁₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: O=C(OCC1=CC=CC=C1)C1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: 323-324°C
Melting pt: 18°C
Solubility: Soluble in Water
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Rokade K, Mali G. World J Microbiol Biotechnol. 2014 Mar;30(3):827-33.

Note: This product is for in vitro research use only