

Cyclophilin F Polyclonal Antibody

Catalog No: #40808

Package Size: #40808-1 50ul #40808-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclophilin F Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Cyclophilin F Polyclonal Antibody detects endogenous levels of Cyclophilin F protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Cyclophilin F.
Target Name	Cyclophilin F
Other Names	PPIF; CYP3; Peptidyl-prolyl cis-trans isomerase F; mitochondrial; PPlase F; Cyclophilin F; Rotamase F
Accession No.	Swiss-Prot: P30405NCBI Gene ID: 10105
Uniprot	P30405
GeneID	10105;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using Cyclophilin F Polyclonal Antibody



Western Blot analysis of K562 cells using Cyclophilin F Polyclonal Antibody

Note: This product is for in vitro research use only