

Triacetin


Catalog No: #S0730

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Triacetin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	218.2
Formulation	C ₉ H ₁₄ O ₆
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CC(=O)OCC(OC(=O)C)OC(=O)C
Pathways: Others
Receptor: Others
Boiling pt: 258-260°C
Melting pt: 3°C
Solubility: Water: 70 mg/mL
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Tsen AR, et al. Int J Y. 2014 Mar 15;134(6):1300-10.

Note: This product is for in vitro research use only