

Delphene


Catalog No: #S0731

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Delphene |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | AChR inhibitor |
| Calculated MW | 191.27 |
| Formulation | C12H17NO |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Nervous system
SMILES: CCN(CC)C(=O)C1=CC(C)=CC=C1
Pathways: Neuroscience
Receptor: AChE
Boiling pt: 111°C
Melting pt: -45°C
Solubility: Soluble in Water
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Xu P, et al. Proc Natl Acad Sci U S A. 2014 Nov 18;111(46):16592-7.

Note: This product is for in vitro research use only