

## Benzyl alcohol

Catalog No: #S0732


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Benzyl alcohol
Brief Description	Inhibitors
Purification	99.34%
Target Name	Dehydrogenase
Calculated MW	108.14
Formulation	C7H8O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:ImmunologySMILES:OCC1=CC=CC=C1Pathways:MetabolismReceptor:ALDHBoling pt:205CMelting pt:-15CSolubility:Soluble in WaterAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Taylor PM, et al. J Feline Med Surg. 2012 Aug;14(8):516-26.

Note: This product is for in vitro research use only