

Pindolol

Catalog No: #S0735


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pindolol
Brief Description	Inhibitors
Purification	99.94%
Target Name	5-HT Receptor agonist; Adrenergic Receptor antagonist
Calculated MW	248.32
Formulation	C ₁₄ H ₂₀ N ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: CC(C)NCC(O)COC1=CC=CC2=C1C=CN2 Pathways: GPCR/G Protein;

Neuroscience Receptor: 5-HT_{1A}; β -adrenergic receptor Boiling pt: Melting pt: 167-171°C Solubility: Soluble in Water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Joseph SS, et al. Naunyn Schmiedeberg Arch Pharmacol. 2003 Dec;368(6):496-503.

Note: This product is for in vitro research use only