

## Sulfamethazine

Catalog No: #S0741


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sulfamethazine
Brief Description	Inhibitors
Purification	99.73%
Target Name	Others
Calculated MW	278.34
Formulation	C <sub>12</sub> H <sub>14</sub> N <sub>4</sub> O <sub>2</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: CC1=CC(C)=NC(NS(=O)(=O)C2=CC=C(N)C=C2)=N1  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: 197°C  
Solubility: DMSO: 56 mg/mL (201.2 mM)  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Fukahori S, et al. J Environ Manage. 2015 Jul 1;157:103-10.

Note: This product is for in vitro research use only